

# CURRY NetStreaming Protocol Quick Guide

This document contains a brief overview of the process of connecting to and retrieving data from a CURRY instance set up as NetStreaming server running NetStreaming protocol version 803.

Data types, command types and functions mentioned in this document refer to the C++ NetStreaming example, which is part of the CURRY installation and by default be found in *C:\Program Files\Neuroscan\Curry 8\Acquisition\CURRY NetStreaming Demo*.

## Setup CURRY as a NetStreaming Server

To set up CURRY as a NetStreaming Server, please refer to chapter *Acquisition > Acquire Parameter Dialogs > Amplifier Control > Configure as NetStreaming Server or Client* in the CURRY User Guide (*Help > User Guide*).

## Use NetStreaming Server Simulator

If no instance of CURRY is available, the *CURRY NetStreaming Server Simulator* tool can be used to simulate the connection process and retrieve a simulated data stream.

To start the server simulator tool, run *C:\Program Files\Neuroscan\Curry 8\Acquisition\CURRY NetStreaming Demo\CURRYNetStreamingServerSimulator.exe* and click **Enable**.

The server simulator operates on port 4455 and provides access to a data stream of 10 channels, containing sine waveforms of different wavelengths.

**Send Event** generates a data packet containing event data.

**Send Impedances** generates a data packet containing impedance data.

The server simulator understands all commands a client can send, but does not react to all of them. For example, commands to start/stop a recording are acknowledged, but no file is being written.

## Commands between Server and Client

Information about commands and the structure of packets can be found in *packets.h*.

## Connecting to and retrieving data from a NetStreaming Server

This section details the process of connecting to a NetStreaming server, requesting information, sending commands and handling data packets.

The connection process is performed in function `CCURRYNetStreamingDlg::OnBnClickedButtonConnect(...)` of the C++ NetStreaming demo project.

General handling of data packets from the server is shown in `CnetStreamingReceiver::HandleMessage(...)`.

- (optional) Start CURRY with parameter '/CS'  
`C:\Program Files\Neuroscan\Curry 8\x64>Curry8.exe /CS`  
This will start CURRY and have it create an unfiled study, so a client can directly connect.
- In client, open socket to server - default port 4455
- Check NetStreaming version of server - optional, but recommended
  - send request `Request_Version` to server
  - wait for reply `InfoType_Version` from server  
body type: long
- Get amplifier status - can be done at any time
  - send `Request_StatusAmp` to server
  - wait for reply `InfoType_StatusAmp` from server  
body type: long
- Set recording file name
  - ensure that file does not exist yet (CURRY will not overwrite a file without user interaction)
  - send request `Request_SetRecPath` to server with filename in body (type `wchar_t[260]`)
  - (server does not send a reply)
  - Note: The recording file name is only accepted for the next acquisition. If CURRY is not connected to an amplifier yet, the requested file patch will be set immediately. If CURRY is currently running an acquisition, the requested file path will be remembered and set by CURRY after the amplifier has been stopped (Stop Amplifier).
- Start amplifier if it is not running yet
  - send request `Request_AmpConnect` to server
  - wait for reply `Server_AcquisitionStart` from server - optional
- Get basic acquisition information from server
  - send request `Request_BasicInfoAcq` to server
  - wait for reply `InfoType_BasicInfoAcq` from server  
body type: `BasicInfoAcq`
- Get channel information from server
  - send requested `Request_ChannelInfo` to server
  - wait for reply `InfoType_ChannelInfo` from server  
body type: one or multiple `NetStreamingChannelInfo`
- Request server to start streaming data
  - send request `Request_StreamingStart` to server
  - if amplifier is running, server will start sending data packets
    - for handling of data blocks containing `DATA_Eeg`,  
see `CnetStreamingReceiver::DataReady(...)`
    - for handling of data blocks containing `DATA_Impedances`,

- see `CnetStreamingReceiver::ImpedanceDataReady(...)`
  - for handling of data blocks containing `DATA_Events`, see `CnetStreamingReceiver::EventDataReady(...)`
- Request server to start impedance check
  - send request `Request_ImpedanceStart` to server
  - server sends `Server_ImpedanceStart` and start sending impedance data packets (`DATA_Impedances`)
- Request server to stop impedance check
  - send request `Request_ImpedanceStop` to server
  - server sends `Server_ImpedanceStop` when impedance check has stopped
- Request server to start recording
  - send `Request_RecordingStart` to server
  - server sends `Server_RecordingStart`
- Request server to stop recording
  - send `Request_RecordingStop` to server
  - server sends `Server_RecordingStop`
- Request server to stop streaming data
  - send `Request_StreamingStop` to server
  - server stops sending EEG, Event and Impedance data packets
- Request server to stop amplifier
  - send request `Request_AmpDisconnect` to server
  - server sends `Server_AcquisitionStop` when amplifier has stopped